

## Thermo Scientific HERAcell Vios 250i Double Stack CO2 Incubators



Copper chamber, double-stack, 255 L / 9.0 cubic feet

Rating: Not Rated Yet

**Price**

Sales price \$10,500.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

Thermo Scientific Heracell VIOS 250i copper chamber co2 incubator, double-stack, standard thermal conductivity, 255 L / 9.0 cubic feet, 120 volts, 50/60 hertz.

With iCAN touchscreen control. TC CO2 sensor.

Relative humidity: =93% at 37°C

Temperature Range (Metric): Ambient +3° to 55°C

Interior: Copper

Shelves: 3 per unit

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**