

**Thermo Scientific Heratherm Oven, 2.2 cu ft.**





Temperature Range (Metric): 50° to 330°C

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Model OMH-60S

2 shelves

Effectively protects samples with an extra layer of security for applications where process reliability and security are paramount. A standard over-temperature alarm and an additional under-temperature alarm ensure that your sample stability. A door lock prevents disruption or tampering of samples, while a door alarm notifies the operator when the door is left open accidentally.

Features

- Broad temperatures range—up to 330°C.
- Advanced digital timer for daily or weekly on / off cycles.
- Stainless-steel interior is easy to clean and corrosion resistant (steel quality AISI 304).
- Programmable controller for temperature ramps and dwells.
- Up to 10 programs can be saved, each with 10 discrete steps.
- Electronically controlled fan speed and damper position can be incorporated.
- Programs can be repeated automatically.
- Access port allows the introduction of sensors for independent data monitoring.
- Simple calibration routine ensures temperature accuracy over time.
- Boost function enables rapid heating —no need to operate the oven 24/7.

Specifications

- Cabinet Material: Coated steel
- Capacity (Metric): 62L
- Dimensions (D x W x H) Interior: 13.5 x 12.9 x 20 in. (34.3 x 32.8 x 50.8cm)
- Dimensions (L x W x H) Exterior: 22.2 x 20.9 x 28.3 in. (56.5 x 53 x 72cm)
- Electrical: Amperage: 12A; Hertz: 60Hz; Voltage: 120V; Wattage: 1440w
- Shelves: 2 supplied/13 max
- Temperature Range (Metric): 50° to 330°C
- Temperature Control: Programmable Microprocessor control with vacuum fluorescent display
- Timer: weekly / real time / hour
- Temperature Stability: ±0.2°C
- Type: Mechanical convection
- Plug Type Nema 5-15

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office

650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)