

## Thermo Scientific NanoDrop 3300 Fluorospectrometer



Broad spectrum fluorescence analysis, using as little as 1  $\mu$ L of sample volume.

Rating: Not Rated Yet

**Price**

Sales price \$8,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

**Description**

With CPU, monitor, keyboard and mouse.

Source Lamp: LED

Min Wavelength (nm): 260

Max Wavelength (nm): 280

Photometric - Absorbance (A): 0

Max Photometric - Absorbance (A): 30

Detector Silicon photodiode

Concentration: 0-1500 ng/ $\mu$ L (dsDNA); 0-45 mg/mL (BSA)

Materials of Construction: 304 Stainless Steel, quartz fiber

Width (cm): 11.5

Depth (cm): 16

Width (in): 4 1/2

Depth (in): 6 19/64

Power (VDC): 6

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**

**226 Miller Ave.**

**South San Francisco, CA 94080**

650-871-5707 office

650-276-7487 fax

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)