

Thermo Scientific Sorvall ST1R Plus-MD Refrigerated Benchtop Centrifuge



Maximum speed: 5,500 rpm with the TX-200 swinging-bucket rotor.

Rating: Not Rated Yet

Price

Sales price \$6,850.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

The Thermo Scientific Sorvall ST1R Plus-MD is a refrigerated benchtop centrifuge designed for routine laboratory and clinical applications. With its swinging-bucket rotor and adapters, it accommodates common 15 mL and 50 mL conical tubes, while providing temperature-controlled centrifugation and convenient programmable operation.

Key Features/Specifications

Maximum speed: 5,500 rpm with the TX-200 swinging-bucket rotor.

Maximum RCF: 5,580 × g with the TX-200 rotor.

Refrigerated: Temperature setting range of approximately –10°C to +40°C.

Swinging-bucket capacity: Up to 4 × 400 mL with the TX-400 rotor; the cell-culture configuration can accommodate up to 16 × 50 mL or 36 × 15 mL tubes.

Convenient controls & safety: Microprocessor control, 9 acceleration / 10 braking profiles, up to 6 stored programs, SMARTSpin imbalance detection, and tool-free Auto-Lock rotor exchange.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com