

## Thermo Scientific TSG30RPGA General Purpose Laboratory Refrigerator



30 cu ft class

Rating: Not Rated Yet

### Price

Sales price \$2,950.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

### Description

Thermo Scientific TSG30RPGA is a TSG Series general-purpose laboratory refrigerator with a single double-pane glass door. It is for everyday lab storage of reagents, media, and other non-critical samples at refrigerator temperatures—not a pharmacy/vaccine or medical-device unit. White-painted cabinet; factory setpoint +5 °C. Model code: TSG series, 30 size class, R refrigerator, P painted, G glass door, A 115 V / 60 Hz.

### Key Features / Specifications

Capacity and size — Rated 27 cu ft / 764 L (sold as the 30 cu ft class). Exterior about 32 × 30.75 × 81 in (D × W × H); interior about 28.5 × 27 × 61.75 in. Shipping weight ~386 lb.

Temperature — Operating range +2 to +10 °C, factory setpoint +5 °C. Digital microprocessor control. Adaptive defrost. Hydrocarbon refrigerant R290.

Glass door and interior — One double-pane glass door, keyed lock, LED lighting. 4 full shelves (about 100 lb / 45 kg each). Rear access port. White painted exterior/interior.

Alarms and quiet operation — High/low temperature, door-ajar, and power-failure alarms (battery backup on the controller). Remote alarm contacts. Whisper-quiet, about 51 dBA. ENERGY STAR and cULus.

Electrical / warranty — 115 V, 60 Hz, NEMA 5-15P, ~3 A. Made for general lab use, not FDA-registered medical storage.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)