

Thermo Scientific TSX5005GA High-Performance Biomedical Lab Refrigerator





51 cu ft., Glass Doors, 120V

Rating: Not Rated Yet

Price

Sales price \$5,250.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

The Thermo Scientific TSX5005GA high-performance biomedical lab refrigerator has a 51 cu ft (1,447 l) capacity, comes with a double glass and runs on 115V power. Designed to safely store medical and laboratory products, including vaccines, pharmaceuticals, media and reagents requiring a temperature range of 2°C to 8°. Whisper quiet with a sound level of 52 dBA, enabling end users to work alongside the refrigerator without distraction. Unique V-drive technology is designed to detect usage patterns such as door openings when a higher compressor speed is needed and periods of stability where the compressor runs at a lower speed, saving energy without compromising protection.. Operating temperature range of 2° to 8°C, factory setpoint at 5°C. Easy to use microprocessor control with interior light, setpoint security and alarms. Standard door locks provide comfort and assurance. Ergonomic door handle.

Specifications:

Capacity: 51 cu ft. (1,447l)

Temperature Settable Range (Setpoint): 3° to 7°C (5°C)

Interior Dimensions (DxWxH): 28.5 x 52 x 58 inches (72.3 x 133.4 x 147.3 cm)

Exterior Dimensions (DxWxH): 37.8 x 56.5 x 78.5 inches (96 x 143.5 x 199.4 cm)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com