

## Thermo Scientific Variomag Biosystem 4 Direct Magnetic Stirrer



Speed Range: 5-120 rpm; Capacity: Up to 2L per position

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

#### Description

The Thermo Scientific Variomag Biosystem 4 Direct is a specialized magnetic stirrer designed for gentle, low-shear agitation of thermally sensitive cell cultures, including microcarrier cultures, culture broths, and cell suspensions, within CO2 incubators. Its hermetically sealed stainless steel housing resists microorganisms and aggressive media, while producing minimal heat to protect sensitive cultures. The unit features a built-in controller with adjustable speeds from 5 to 120 rpm, offering smooth acceleration and reliable performance for up to 2L per position across four stirring positions.

#### Specifications:

Model: Variomag Biosystem 4 Direct

Stirring Positions: 4

Speed Range: 5-120 rpm

Capacity: Up to 2L per position

Material: Hermetically sealed stainless steel housing

Voltage: 100-240V

Power: 1.5W per position

Features: Low-shear agitation, minimal heat output, built-in controller

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)