

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Thermo Sorvall Legend XTR Refrigerated Benchtop Centrifuge



15,200 rpm maximum speed, 25,314xg maximum g-force

Rating: Not Rated Yet

Price

Sales price \$7,500.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Features simple operation with intuitive control panel and digital display, auto-locking rotor system with push-button installation and rotor exchange, advanced rotor management system. Energy efficient design.

Catalog no. 75004521

15,200 rpm maximum speed, 25,314xg maximum g-force

4x750mL maximum capacity

-10 to +40°C temperature range

120V, 60Hz, 12A

With Thermo TX-1000 swinging bucket rotor Catalog no. 75003017 and 4 buckets, two 50mL tube adapters, two 15mL tube adapters, and four biocontainment lids, 4,200 rpm maximum speed with TX-1000 rotor.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com