

## Thermo Titerplate Shaker



## Shakers/Stirrers/Hotplates: Thermo Titerplate Shaker

---

Speed Range: 40-1100 rpm

Rating: Not Rated Yet

**Price**

Sales price \$550.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

The Thermo Titerplate Shaker Model 4625 is a reliable benchtop shaker designed for efficient mixing of up to four 96-well titer plates, ideal for virology, immunology, serology, and hematology applications. It features a gentle orbital motion with a 0.12-inch (3.2mm) orbit and solid-state speed control, allowing adjustable speeds from 40 to 1100 rpm. The built-in timer supports continuous operation or timed runs from 30 seconds to 5 minutes, with an automatic shut-off at the end of timed cycles, ensuring consistent and repeatable results.

### Specifications:

Model: 4625

Plate Capacity: Up to 4 standard 96-well titer plates

Speed Range: 40-1100 rpm

Orbit Size: 0.12 inches (3.2mm)

Timer: 30 seconds to 5 minutes (adjustable in 30-second increments) or continuous

Max Load: 2 lbs (0.9 kg)

Operating Temperature: 5°C to 40°C (41°F to 104°F)

Features: Automatic shut-off, solid-state speed control, stainless steel retaining springs

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)