

Torrey Pines EchoTherm Chilling Incubator



The EchoTherm™ IN30 Series bench-top, solid-state convection chilling/heating incubators is reliable, accurate and easy to use. Peltier-based for both heating and chilling, the units have no compressors or CFCs, making them truly energy efficient.

Rating: Not Rated Yet

Price

Sales price \$1,150.00

[Ask a question about this product](#)

Manufacturer [Torrey Pines Scientific](#)

Description

• **!!!! - RECENTLY SOLD - !!!!**

-
- Peltier-based. No compressors, no CFCs. Energy-efficient
- Temperature range from 4.0°C to 70.0°C
- Fully programmable or non-programmable controls
- 27.5 L volumes
- Timer with alarm and Auto-Off
- Electronic calibration and power failure protection
- RS232 interface
- 12 month warranty
- UL, CSA and CE listed

Applications Include:

- Protein crystal growth
- BODs
- Culture growth above, below or at ambient
- Enzyme reactions and deactivations
- Incubating marine cultures below ambient
- Ligations at 14.0°C–16.0°C
- Hybridizations
- Storing Oocytes at 17.0°C
- Storing DNA libraries