

## Turner BioSystems Modulus Fluorometer



SOLD A multifunctional, single-tube instrument for fluorescence, luminescence and absorbance.

Rating: Not Rated Yet

**Price**

Sales price \$2,650.00

[Ask a question about this product](#)

Manufacturer [Turner BioSystems](#)

### Description

**SOLD**

Designed to give you the utmost flexibility for measuring fluorescence, luminescence and absorbance. Modulus allows you to buy what you need now, and to add what you want later.

- Operates as a dedicated fluorometer, or a combination fluorometer/luminometer
- Optional Absorbance Module permits absorbance measurements in the visible wavelength spectrum

Available with Green or Blue filter modules.