

Vacuubrand BVC Professional Fluid Aspiration System



Multiples in stock.

Rating: Not Rated Yet

Price

Sales price \$1,375.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

Description

Features:

Powerful and long-lasting chemistry diaphragm pump for universal economical use
Extremely compact, space saving and low noise and thus the ideal tool in the workplace
Easy operation due to electronically adjustable suction power via touch panel, for strong or sensitive suction
Autoclavable collection bottle designed to reduce aerosol and foam formation; with hydrophobic protection filter 0.2 μm , autoclavable quick-couplings available as accessory
Ready to accept a second aspiration hand controller (VHCpro) to support two users

Specifications:

Vacuum Pumps: Vacuubrand BVC Professional Fluid Aspiration System

Max. pumping speed 50 Hz: 0.7 m³/h (Vacuum pump)
Max. pumping speed 60 Hz: 0.5 cfm (Vacuum pump)
Ultimate vacuum: 150 mbar / 112 torr
Vacuum controller: Touchpanel
Control range: 150 - 850 mbar (adjustable)
Collection bottle: 4l Polypropylene
Max. back pressure (abs.): 1.1 bar
Inlet connection: Suction controller
Outlet connection: Hose nozzle DN 8 mm / silencer
Rated motor power: 0.04 kW
Rated motor speed at 50/60 Hz: 1500/1800 min⁻¹
Protection class: IEC 60529IP 40
Weight: 7.7 kg
Noise (sound pressure level) in dBA at 50 Hz or 1500 rpm/62% (VARIO)/1500 rpm (VARIO-SP)/12500 Upm (VACUU·PURE®): 49 dBA
NRTL certification: No
Rated mains voltage range 1: 15 V
Mains frequency 1: 50-60 Hz

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com