

## Vacuubrand PC 3001 Vario Select - Appears Unused



Ultimate vacuum 2.0 mbar (1.5 torr) / 4 mbar with gas ballast; max pumping speed 2.0 m<sup>3</sup>/h (50 Hz) / 1.18 cfm (60 Hz)

Rating: Not Rated Yet

**Price**

Sales price \$3,950.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

### Description

The Vacuubrand PC 3001 VARIO select is a compact, oil-free chemistry vacuum pumping unit with integrated VACUU-SELECT controller. It combines a speed-controlled 4-head/3-stage diaphragm pump with intelligent vacuum control for precise, automated process optimization—ideal for rotary evaporation, vacuum drying, concentration, and other lab applications involving even high-boiling solvents. It features excellent solvent recovery via a glass inlet separator and exhaust condenser.

### Key Features

Intelligent VARIO Speed Control + VACUU-SELECT: Automatic boiling point detection, continuous vacuum optimization, and application-based interface (predefined methods + drag-and-drop custom apps) for walk-away operation and superior process control.

High Performance & Chemical Resistance: Ultimate vacuum 2.0 mbar (1.5 torr) / 4 mbar with gas ballast; max pumping speed 2.0 m<sup>3</sup>/h (50 Hz) / 1.18 cfm (60 Hz); 4 heads / 3 stages for reliable handling of aggressive vapors.

Quiet & Efficient: Extremely low noise (42 dBA), minimal vibration, and low power consumption (0.16 kW motor); long maintenance intervals thanks to VARIO technology.

## Vacuum Pumps: Vacuubrand PC 3001 Vario Select - Appears Unused

---

Excellent Solvent Recovery: Includes protected glass inlet separator (AK) and exhaust vapor condenser (EK) for nearly complete solvent capture using chilled water or recirculator.

Compact & User-Friendly: Small footprint (303 × 306 × 400 mm, 8.2 kg), intuitive touchscreen-like controller, multiple interfaces (Ethernet, USB, RS-232), and easy setup for rotary evaporators and other systems. (Note: Oil-free, chemically resistant design with ATEX conformity for internal explosive atmospheres.)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)