

## Virtis VirSonic 100 Cell Disruptor



With fine-tip probe.

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

### Description

**Generator:**

Dimensions: 33cm W x 19cm D x 17cm H (13" x 7.5" x 6.7")

Weight: 4.2 kg (9.2 lb.)

Input Voltage: 115 VAC nom. (switchable to 220 VAC nom.) 50/60 hz

Full Load Current: 1.65 amps at 115v/ 0.8 amps at 220v

Fuse Rating and Type: 3 amp/1.5 amp Fast acting

Voltage Tolerance: +/- 10% nominal voltage

Output Voltage/Frequency: 950 Vrms (maximum) 22.5 kc (nom)

**Convertor:**

Dimensions: 17cm L x 3cm Dia. (6.7" L. x 1.18" Dia.) without probes Weight: .2 kg (.5 lb.)

Materials: Aluminum case and Front Driver

**Environmental:**

Pollution Degree 1

Temperature Limits: 50° F - 110° F (10° C - 43° C)

Barometric Pressure: Unlimited (Note: In high vacuum areas, additional cooling provisions may be needed.)

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**sales@biotechequipmentsales.com**

