

## VWR 1130S Heating Circulating Water Bath 6L



Appears Unused: 120V / 60Hz / 11 A; Temp range: +5-150 C

Rating: Not Rated Yet

**Price**

Sales price \$750.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

This six-liter capacity heating circulating bath appears unused. It has a 2-speed Simplex Pump with a pressure flow rate of 7 or 15 liters/min.

Electrical req's: 120 V / 60 Hz / 11 Amps;

Temp range: +5-150 C

Dimensions: (L x W x D): 14-3/4 x 8-1/4 x 14 in. (37.5 x 21 x 35.5 cm)

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)