

VWR 1160S Chiller



Ideal for direct immersion of samples or circulating through external systems. The standard controller features a two-speed pressure pump and an easy-to-set temperature encoder with three user-set presets to allow for rapid setpoint change. The analog con

Rating: Not Rated Yet

Price

Sales price \$1,200.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

!!! SOLD !!!

The VWR 1160S refrigerated/heating circulating bath has a reservoir capacity of 6 liters. Per its technical specifications, it can achieve temperature, ranging from -20°C to 150°C, with stability of $\pm 0.05^\circ\text{C}$.

This unit has a stainless steel interior and lid.

Please call us at (650) 871-5707 to check availability of this product!