

VWR 1680 High-Performance Oven



NEW VWR 1680 High-Performance Mechanical Oven, 14.6 cu/ft, double door, 40 to 260°C

Rating: Not Rated Yet

Price

Sales price \$8,500.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

!!!! - RECENTLY SOLD - !!!!

This VWR 1680 High-Performance Mechanical Oven is BRAND NEW in its original packaging.

It has two doors, 6 shelves, the interior is 42" W x 20" deep x 30" high, has a Programmable 24-Step Ramp and Soak Controller, a True Horizontal Airflow Design, Type-304 Stainless Steel Interior and Exterior, and Three-inch Intake and Exhaust Ports.

Temperature range: 40 to 260°C 208–230V, 50/60Hz, 6300W, 26.5A