

VWR 1927



VWR 1927 40 cu/ft CO2 incubator w/Wheaton 11-tiered roller rack, NEW

Rating: Not Rated Yet

Price

Sales price \$8,650.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

!!! SOLD !!!

This VWR 1927 40 cu/ft CO2 incubator w/Wheaton 11-tiered roller rack has microprocessor temperature control, IR CO2 control, 4 interior electrical outlets, all stainless interior, and has a flip-out ramp for the roller apparatus.

The Wheaton R2P 11-place roller apparatus is microprocessor controlled, hold up to 55 roller bottles, accurately displays bottle rotation speed and programmable soft start, stop, and reverse modes. 0.25 to 8.75 rpm bottle speed.

Recently Sold: VWR 1927

Both are BRAND NEW!!!

Please call us at (650) 871-5707 to check availability of this product!