

VWR 5706



This -80 upright ultralow freezer is 23 cu/ft, 230V and manufactured by ThermoForma & is 6 years old. Current model.

Rating: Not Rated Yet

Price

Sales price \$5,500.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

!!! SOLD !!!

This -80 upright ultralow freezer is 23 cu/ft, 230V, manufactured by ThermoForma & is 6 years old. It is the current model.

Its major features are:

- Microprocessor Temperature Controller
- Audible/Visual Alarm System
- Automatic Voltage Compensator
- Cascade Refrigeration System
- Non-CFC, Non-HCFC, Nonflammable Refrigerants

Providing reliable performance, temperature stability, and long-term storage protection, these freezers are ideal for industrial and biomedical storage applications where sophisticated controls and alarm systems are not required. Freezers offer an affordable answer to long-term storage needs at stable temperatures.

These freezers feature a microprocessor temperature control system. Air-cooled, hermetically-sealed dual compressors work in unison and cycle on demand of the microprocessor temperature controller (-40°C units have a single-stage compressor). The automatic voltage compensator detects line frequency, monitors incoming voltage, and adjusts compressor voltage to reduce wear and tear on the compressors. Front-to-back air circulation system results in cooler compressor temperatures, increased performance, and longer compressor life. Easy-to-read, green LED continuously displays chamber temperature. The probe for the temperature controller is located inside the freezer chamber to show actual temperature of stored products.

Audio/visual alarm system alerts operator to over-temperature, under-temperature, power failure, door ajar, and low battery conditions. Battery back-up provides power in the event of a power failure.

Modern design includes an easy-to-grasp handle for one-handed operation. Dual wheel, locking swivel casters allow mobility as well as secure positioning. Hinged grille swings out, exposing the foam air filter for easy cleaning or replacement. Industrial-grade cabinets are constructed of high-quality, heavy-gauge, cold-rolled steel. The powder-coated finish is rust- and chip-resistant.

Please call us at (650) 871-5707 to check availability of this product!