

## VWR 9000L Relative Humidity Chamber

PRODUCT  
Image  
Unavailable

PRODUCT  
Image  
Unavailable

VWR 9000L Relative Humidity Chamber, refrigerating, 5.5 cu/ft, 10C to 70C, 10% above ambient to 95% RH

Rating: Not Rated Yet

**Price**

Sales price \$6,250.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

**!!! SOLD !!!**

VWR 9000L Relative Humidity Chamber, refrigerating, 5.5 cu/ft, 10C to 70C, 10% above ambient to 95% RH, has LED readouts in temp & RH, over-temp protection, 19-gallon water jacket, stainless steel interior & shelves, very clean.

**Please call us at (650) 871-5707 to check availability of this product!**