

## VWR Analog Rocker - 1 Tier



Platform Size: 14 x 11 inches

Rating: Not Rated Yet

**Price**

Sales price \$525.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

**Description**

The VWR Analog Rocker is designed for gentle mixing of samples, providing a reliable and compact platform for various applications. With control knobs for speed, tilt, and time, the user can easily make adjustments while the unit is in motion.

**Specifications:**

Platform Size: 14 x 11 inches (35.6 x 27.9 cm)

Speed Range: 0 to 60 RPM

Orbital Motion: 1 inch (2.5 cm)

Power: 115V, 60Hz

Applications: Perfect for mixing solutions, staining and destaining gels, hybridization procedures, hematology, suspending cells, and gentle agitation in biochemistry and molecular biology tasks.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)