

VWR AS12



VWR AS12 benchtop microprocessor controlled autoclave. Sterilizes up to three 2L flasks in one cycle.

Rating: Not Rated Yet

Price

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

SOLD

The VWR AS12's design is microprocessor based, permitting one touch operation. It features four preset programs for liquids, wrapped and unwrapped goods, and packs. Also, custom programs can be entered for liquid or hard goods. The LCD display shows the program, run status, and any nonconformity. The built in water pump assures accurate steam for every cycle. Can fit three 2L flasks. Other safety features include door interlock so that the door cannot be opened until chamber depressurization, and an access port to permit verification of temperature and pressure. The external printer (shown here) records sterilization performance when required.