


VWR Cryopro L-30 Liquid Nitrogen Tank -- Brand New!



Overall Height x Outside Diameter, cm (in.): 61.1 x 43.2 (24 x 17)

Rating: Not Rated Yet

Price

Sales price \$675.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

Specifications for VWR CryoPro Liquid Dewars, L Series L-30 Liquid Dewar:

Description: L-30 Liquid Dewar

Static Evaporation Rate, L/Day: 0.22

Liquid Nitrogen Capacity, L: 32

Weight Empty, kg (lbs.): 12 (26.4)

Weight Full, kg (lbs.): 38 (83.8)

Neck Diameter, cm (in.): 6.4 (2 1/2)

Overall Height x Outside Diameter, cm (in.): 61.1 x 43.2 (24 x 17)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com