

## VWR Digital Heat Block



These units are ideal for incubation and activation of cultures, enzyme reactions, and a wide variety of other laboratory procedures.

Rating: Not Rated Yet

### Price

Sales price \$200.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

**!!! SOLD !!!**

These units are ideal for incubation and activation of cultures, enzyme reactions, and a wide variety of other laboratory procedures. Front-mounted controls with independent LED displays show actual and set point temperatures. Microprocessor offers a stability of 0.1°C, and independent controls allow adjustments of the temperature settings in  $\pm 0.1^\circ\text{C}$  increments. A timer will display elapsed time or when programmed to user-defined limit, will count down to zero and shut off unit. Audible alarm will sound when time reaches zero.

**Please call us at (650) 871-5707 to check availability of this product!**

