


VWR Digital Rocker



Speed range: 1 to 50 rpm

Rating: Not Rated Yet

Price

Shakers/Stirrers/Hotplates: VWR Digital Rocker

Sales price \$675.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

Cat # 12620-906

12.75 x 10" platform

Speed range: 1 to 50 rpm

Speed accuracy: 1 rpm

Tilt angle: 1 to 5

Timer: 1 sec to 9,999 min

Maximum load: 4.5 kg

Overall dimensions (W x H x D): 11 x 5 x 17 in

Electrical requirements: 120V; 50/60Hz

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com