

VWR Forced Air Oven 3,65CF



Temperature Range: 50C to 250°C

Rating: Not Rated Yet

Price

Sales price \$1,150.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

Forced air convection oven.

Specifications:

Exterior Dimensions: 25.2W x 36.2H x 29.1D"

Interior Dimensions: 18.3W x 27.9H x 21.4D"

Heat Up Time: 25 min.

Recovery Time: 4 min.

Temperature Uniformity: $\pm 3.5^{\circ}\text{C}$ (based on 27 sensors)

Fluctuations: $\pm 0.4^{\circ}\text{C}$

Electrical: 120V

Temperature Range: 50C to 250°C

Temperature Setting Accuracy: 0.1°C

Catalogue #: 89511-412

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com