

VWR Incubating Microplate Shaker



Temperature Range: Ambient +5...65 °C (149 °F)

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

Ideal for cell cultures, hybridizations, and suspensions.

Specifications

Temperature Range: Ambient +5...65 °C (149 °F)

Temperature Uniformity: $\pm 0.5^{\circ}\text{C}$ at 37°C

Speed accuracy: $\pm 2\%$

Max. Load: 2.3 kg (5 lbs.)

Overall Dimensions: 28W x 19.7H x 43.2D cm (11 x 73/4 x 17")

Timer: 1 second to 160 hours

Shipping Weight: 13.6 kg (30 lbs.)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com