

## VWR Micro Plate Shaker, New in Box!



Low-profile design with a 3 mm orbit

Rating: Not Rated Yet

**Price**

Sales price \$750.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

The VWR Micro Plate Shaker (Cat. # 12620-926) is a compact, digital orbital shaker designed specifically for mixing 96-well and other microplates. It is new in box (as shown), features a low-profile design with a 3 mm orbit, and is ideal for immunoassays, ELISA, cell culture, PCR prep, and general microplate shaking applications. It holds 2 or 4 standard microplates (or deep-well plates) securely without additional clamps.

### Key Features & Specifications:

Speed Range & Accuracy — 100–1,200 RPM with  $\pm 2\%$  accuracy (up to 999 RPM) for consistent, reproducible mixing.

Orbit — 3 mm orbital diameter, optimized for gentle yet effective microplate shaking.

Digital Controls & Timer — Microprocessor-controlled with independent LED displays for speed and time; timer from 1 second to 160 hours (or continuous mode) with audible end-of-cycle alarm.

Capacity — Holds 2 or 4 microplates (or 2 micro-tube racks); aluminum platform with secure positioning.

Compact & Practical Design — Dimensions: 11" W x 4" H x 17" D (27.9 x 10.2 x 43.2 cm); 120V, 50/60 Hz operation; lightweight (~17–25 lbs shipping weight); quiet operation with over-speed protection.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc. has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)