

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

VWR model VWR-HC-FRP-20 - 20 cu. ft. Digital Laboratory Flammable Storage Refrigerator from 2023



Temperature range: Adjustable 1°C to 10°C

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

The VWR-HC-FRP-20 is a 20 cu. ft. upright laboratory refrigerator designed for safe cold storage of flammable and combustible materials (in sealed containers) in unclassified general laboratory work areas. It is a basic/standard-series digital model using natural hydrocarbon refrigerant, compliant with relevant NFPA guidelines.

Key Features/Specifications

Capacity: 20 cu. ft. (?566 L)

Temperature range: Adjustable 1°C to 10°C, managed by a basic digital microprocessor controller with digital temperature display (°C/°F switchable) and sample-simulated readout (0.1° resolution on many equivalents)

Alarms & monitoring: Audible/visual high/low temperature alarms; remote alarm contacts; often includes battery backup support for the display/alarm module and dual probes (air + sample bottle)

Door & access: One solid right-hinged swing door (not reversible), keyed lock, magnetic gasket, heavy-duty handle; probe access port (typically ?? or ¾?)

Interior: 5 adjustable shelves (usually 4 full + 1 half); forced-draft/forced-air circulation for uniformity and recovery

Refrigerant & defrost: Eco-friendly natural hydrocarbon (R600a); cycle defrost

Construction & safety: White textured/powder-coated steel exterior; high-density urethane foam insulation; utilizes non-sparking and/or ATEX-rated components for an intrinsically safer interior; meets NFPA guidelines 45 and 70 (or 30) for flammable materials storage in unclassified locations; ETL/C-ETL listed (to UL 471 standards)

Dimensions (exterior, approx.): 30.75" W × 29" D × 74.625" H (interior roughly 24.75" W × 20.75" D × 65.25" H)

Electrical: 115 V, 60 Hz, 22.5–3 A, ¼ HP; NEMA 5-15 plug

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com