

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

VWR Vertical SpaceSaver Glassware Washers (Doublestack)



New. Top DI rinse/standard rack and bottom spindle rack ready.

Rating: Not Rated Yet

Price

Sales price \$12,995.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

Maximum capacity and flexibility in a small footprint for laboratory glassware cleaning. These vertical SpaceSaver glassware washers deliver high-volume, reliable washing and drying of lab glassware in a compact design ideal for space-constrained laboratories. Independent microprocessor controls, flexible rack options, and efficient performance while using significantly less floor space than comparable-capacity units.

Key Features:

- Up to 33% smaller footprint than washers of similar capacity (vertical/stackable design).
- Single- or double-washing chamber models with up to four loading levels total.
- Independently microprocessor-controlled chambers for maximum throughput and flexibility.
- Premium 304 stainless steel construction for durability and corrosion resistance.
- Soft-touch controller with multiple cycle selections (e.g., rinse only, normal, daily, heavy wash, timed/custom options).
- Temperature selection up to 70 °C (158 °F) for optimal cleaning of different loads.
- Enhanced dual-fan drying system that helps prevent impurities and contamination.
- Noise-buffering insulation for quieter operation.
- High-efficiency water management system for energy and water conservation.
- Spindle-rack-ready options (on applicable models) plus included top/bottom stainless-steel open racks; compatible with a wide range of accessories; supplied with mounting feet and wheels.

Miscellaneous: VWR Vertical SpaceSaver Glassware Washers (Doublestack)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com