

Watson-Marlow, Inc PumpPro DPM



Microprocessor control delivers rapid, precise dispensing or transfer of up to three litres per minute with a dispensing accuracy better than +/- 1%.

Rating: Not Rated Yet

Price

Sales price \$900.00

[Ask a question about this product](#)

Manufacturer [Watson-Marlow, Inc](#)

Description

Microprocessor control delivers rapid, precise dispensing or transfer of up to three litres per minute with a dispensing accuracy better than +/- 1%. Membrane keypad enables operator to easily specify required volume, number of dispenses, and pump speed. Select digital display as either flowrate (mL/sec.) or speed (rpm). Dispensing may be initiated by specified time interval or with a foot switch or hand switch. Manual control available for continuous metering or transfer. Adjustable "soft-start" can be programmed to reduce splashing; end-of-dose reversing capability prevents dripping. Includes standard size, easy loading pumphead with lockable guard. Features simple calibration as well as recalibration without batch interruption.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com