

Welch DryFast



Welch DryFast Model 2025 self cleaning vacuum system perfect for rotary evaporators

Rating: Not Rated Yet

Price

Sales price \$1,250.00

[Ask a question about this product](#)

Manufacturer [Welch](#)

Description

!!! SOLD !!!

This Welch DryFast Model 2025 self cleaning vacuum system is perfect for rotary evaporators. It is resistant to chemical attack, even from the most aggressive vapors and containing fluorinated plastics on all wetted surfaces.

Residue accumulation on the diaphragms & valves is reduced by a glass inlet separator that inhibits the ingesting of liquid & particulate contaminants which can impede performance. A gas ballast/vent switch acts to minimize vapor condensation inside the pump, and an automatic self-cleaning purge rids the pump of unwanted residue after every run. Very clean.

Please call us at (650) 871-5707 to check availability of this product!