


Microscopes: Zeiss Primovert Inverted Phase Contrast Microscope

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Zeiss Primovert Inverted Phase Contrast Microscope



with two WF10x/20 eyepieces, 4-slot objective nosepiece turret equipped with x4 Zeiss Primo objectives

Rating: Not Rated Yet

Price

Sales price \$3,995.00

[Ask a question about this product](#)

Manufacturer [Zeiss](#)

Description

Ideal for routine cell and tissue culture observations, system includes a binocular headpiece with two WF10x/20 eyepieces, 4-slot objective nosepiece turret equipped with x4 Zeiss Primo objectives:

- 1) Plan-ACHROMAT 4x/0.10
- 2) Plan-ACHROMAT 10x/0.25 Ph1
- 3) LD Plan-ACHROMAT 20x/0.30 Ph2
- 4) LD Plan-ACHROMAT 40x/0.50 Ph2

LD condenser 0.4 (working distance 55mm), 2-position phase contrast slider with Ph1/0.4 and Ph2/0.4 apertures, 6V 30W halogen illumination, 12VDC.

Comes with our 90-day warranty.

Microscopes: Zeiss Primovert Inverted Phase Contrast Microscope

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com